



# Product Reference Guide


## PNEUPAC™ paraPAC plus™ P310 VENTILATOR

### Quick Reference Guide

## Connect and Set up the Ventilator

### 1. How do I connect the ventilator? (Image A)



- › Connect 1 to wall or cylinder, oxygen only
- › Connect 2 to pressure monitoring line
- › Connect 3 to patient breathing circuit



1. Oxygen supply gas inlet connection
2. Pressure monitoring connection manometer
3. Patient outlet connection

### 2. How do I set up the breathing system and the HME filter? (Image B)


- › Attach patient breathing circuit to the connection (3)
- › Attach the monitor line to the pressure monitoring connection (2)



## Set the controls on the pneupac paraPAC Plus P310 ventilator

### 3. What controls does the ventilator have? (Image C)

- › Set mode to **ON** and select **Position I** for continuous mandatory ventilation (3)
- › Oxygen concentration: select 50% or 100% (7)
- › Frequency knob controls breaths per minute (5)
- › Tidal volume knob alters the gas volume to patient (4)
- › Pressure is displayed on the pressure dial (2)
- › Pressure can be adjusted using the relief alarm pressure control (1)
- › Positive expiratory pressure can be adjusted using the PEEP control (6)
- › The ventilator maintains a 1:2 I:E ratio; 1-part inhalation, 2-parts exhalation. The time of each phase is determined by the flowrate required to deliver a selected number of breaths and tidal volume



**icumedical**  
human connections



# Product Reference Guide

## PNEUPAC™ PARAPAC PLUS™ P310 VENTILATOR

### Quick Reference Guide

## Operating the pneupac paraPAC Plus P310 ventilator

### 4. How do I select the mode?


- Set mode to **ON** and select **Position I** for continuous mandatory ventilation

### 1. How do I set the parameters? (Image C)

- Set **frequency** and **tidal volume** (4 and 5)
- Set parameters appropriate for patient (**yellow** for adults)
- Set **maximum pressure** limit between 20 and 60 cmH<sub>2</sub>O using the relief pressure control (1)
- Positive and expiratory pressure can be adjusted using the **PEEP** control (6)
- Set oxygen percentage to either 50% FiO<sub>2</sub> or 100% FiO<sub>2</sub>, using the Air-Mix control (7)

### 6. How do I check the device?

- Attach the pressure monitor line to the ventilator via the pressure monitoring connection (Image A)
- Once the mode is set, check alarm indicators flash in sequence (Image D), clockwise from high pressure
- Once frequency and tidal volume are set, cover connection port, check peak inflation pressure. The pneumatic alarm sounds and high-pressure alarm shows red



### 7. What do the alarms mean? (Image D)

- Pressure jumps excessively at beginning of inspiration (1)**
  - Airway obstruction
  - Flow rate is too high; decrease tidal volume setting
- Pressure jumps excessively at the end of inspiration (1)**
  - Tidal volume setting is too high; decrease tidal volume or frequency setting or
  - Could be caused by mucous plugging; patient requires suctioning
- Sudden increase in patient airway pressure (1)**
  - Tube kinking, patient coughing, lessening of sedation causing asynchrony or other obstructions
- Inflation pressure too low (3)**
  - Leak in the circuit or disconnection from patient; check hose or connection leaks
  - Check that ventilation parameters



Rx Only. For safe and proper use, refer to the appropriate manual, guide, or instructions for use. Figures are rendered as graphic representations to approximate actual product.